

You will get the most effective results, taking proper care, having balanced meals and using DeVita Cosmo.

WELLNESS AND BEAUTY

HOW IT WORKS

DEVITA COSMO



It combines the effects of bioresonance therapy and quantum radiation technologies. This double action provides the best results and reduced time of use. Optionally, you can use the device without quantum nozzle.



It fights aging signs, lifts the skin and makes it elastic. It contributes to collagen production, improves penetration of cosmetic products into deeper layers of the skin and accelerates regeneration processes.



It eliminates inflammation and redness of the skin by improving circulatory, excretory and hormonal systems' function, as well as function of the gastrointestinal tract. Effective against acne, rosacea and other skin inflammations.



It has a beneficial effect on hair and nails. DeVita Cosmo programs strengthen nail plates and hair, and harmonize the state of all internals, which is positively reflected in your appearance.



UNIVERSALITY
a device which improves the health of the entire body



EFFECTIVENESS
the technologies are recognized by the leading scientific organizations



SIMPLICITY
the interface and the scheme of application are user friendly



SAFETY
contraindications are minimum



PORTABILITY
a comfortable device is always at hand



URGENCY
health, beauty and youth support




WWW.DEVITA.STORE



DeVita-technologies are the fastest, easiest, most reliable and modern way to your health improvement.

ACCORDING TO INTERNET STATISTICS



RECOMMENDATIONS FOR USE

DeVita Cosmo is a cosmetic and revitalizing device that combines bioresonance technologies and quantum properties of colour. Unlike cosmetic products, DeVita Cosmo influences the deeper skin layers and improves it on the inside. The health-improving programs of the device are suitable for the whole family.


DEVITA COSMO IS USED TO:



The unit contains 20 basic programs providing complex effect.

MODE OF ACTION

The DeVita Cosmo action is based on the results of research in the area of quantum light theory. Visible light has a wide colour spectral composition consisting of radiation with different wavelength and frequency ranges. Each colour may have particular effects on certain organs or systems through receptors in the skin. Zakharyin-Ged areas located in various parts of the body, including face and neck, are influenced under the principle of reflex response. Quantum radiation restores the balance of body systems and contributes to their functional restoration due to the reflex action at the nerve endings and to the activation of biologically active substances synthesis.



DeVita Cosmo – is a household wellness device.

PROGRAMS OF DEVITA COSMO

Skin without wrinkles	1	11	Thyroid gland
Facial tonus	2	12	Heart and vessels
Women's health	3	13	Spleen functioning
Men's health	4	14	Key of life
Liver functioning	5	15	Clear eyes
Stressless life	6	16	Clean skin
Stomach control	7	17	Beauty plait
Digestion	8	18	Master of manicure
Like clockwork	9	19	Face without defects
Kidneys stimulation	10	20	Light feet